

 <p><b>VERSATILT OVERBED TABLE</b> 15x30" WALNUT MELANINE TOP 26.5-39" INFINITE HT ADJUST 950H, \$249.95</p>	 <p><b>OVERBED TABLE WITH LOCKING WHEEL STOPS</b> NON-TILT TOP:15.5Wx30L LVA1022, \$149.95</p>	 <p><b>FOOT STOOL, 300LB CAP</b> GUN METAL FINISH 14.25Wx11.25Dx9.25H 6065, \$44.50</p>	 <p><b>FOOT STOOL WITH HANDLE</b> SILVER VEIN FINISH 300LB CAP 6066, \$65.50</p>
---	---	---	---

**PATIENT LIFTS**



- Silver Vein steel construction
- 6-point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base
- 450 LB weight capacity
- Lifetime warranty (frame), 1 year limited (pump)  
13023SV, \$895.00

The PL400HE is Bestcare's answer to home care needs. Featuring a 400 LB maximum weight capacity and Bestcare battery powered control system for power lifting, the PL400HE is easily adaptable for a wide range of challenges while still maintaining a high standard of comfort for patients. The use of Bestcare's battery powered electronic system provides a low-cost solution for power lifting while providing convenience for caregivers. Designed for home care use and effective at providing quality care.



- 400 LB carrying capacity
- Six-point spreader bar
- Electronic control system
- Optional digital weigh scale
- Optional two-point spreader bar  
PL400HE, \$1,975.00

**SLINGS**



MESH SLING FULL BODY		WITH COMMODE OPEN
SL-R110	MED	WEIGHT RANGE 100-195 LBS 450LB CAP
SL-R111	LG	WEIGHT RANGE 150-300 LBS 450LB CAP
WITH COMMODE OPENING:		
SL-R114	MED	WEIGHT RANGE 100-195 LBS 450LB CAP
SL-R115	LG	WEIGHT RANGE 150-300 LBS 450LB CAP
SL-R116	XLG	WEIGHT RANGE 210-400 LBS 450LB CAP

SLING DIVIDED LEG MESH		HOYER TYPE
SL-UM882	MED	FITS WEIGHT RANGE 140-200 LBS
SL-UM883	LG	FITS WEIGHT RANGE 200-450 LBS
HOYER TYPE:		
SL-HC70002	MED	U-SLING 99-210 LBS 500LB CAP
SL-HC70001	LG	U-SLING 150-300 LBS 500LB CAP
SL-HC70000	XLG	U-SLING 270-600 LBS 600LB CAP
SL-HC70012	MED	W/HEAD SUP'RT 99-210 500LB CAP



# FREEDOM

TRANSFER PATIENT LIFT



## ERGONOMIC HANDLES

Comfortable and secure hand grips in the front and rear for easy transport.

01

## SEAT & FRAME LOCK

Dual locking to ensure the safety of the patient and the caregiver.

02

## DUAL CONTROLLED

Features built-in control and a detachable remote for lifting and lowering the patient.

03

## ADJUSTABLE

The seat height may be adjusted from 15.5" to 23.5" to accommodate height differences when transporting.

04

## FOOTREST

Keep the patients feet comfortable and secure when transporting patients.

05



THE  
BED

THE  
TOILET

THE DINING  
ROOM

THE LIVING  
ROOM

YOUR  
FAMILY

**\$2,799.00**

Overall Height: 34.5" – 42.5"  
 Seat Width: 19.5"  
 Seat Height: 15.5" – 23.5"  
 Commode Pail: 6 Qt  
 Required Clearance: 4.5"  
 Rear Wheel Size: 3"  
 Wheel Base Length: 25"  
 Wireless Lithium Battery  
 Footrest Size: 13.75" x 6"  
 Wheel Base Width: 22.3"

**VERSA-HELPER TRAPEZE**



- Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance
  - Rubber bumper guards protect walls
  - Hand bar location fully adjustable
  - Height adjustable from 43" to 52"
  - VersaGuard coated clamps and hand bar
  - For use with metal frame headboards only
  - Assisted weight capacity of 250lbs
- MDS80615T (TRAPEZE ONLY), \$195.95  
 MDS80616T, (FLOOR BASE), \$169.95



Floor Base

**ELECTRIC HOSPITAL BEDS**



- Constructed from durable steel with a spring deck and split-frame design
- Spring deck is adjustable from 15.7"-22.8"
- Motor is compatible with a 9-volt battery to adjust bed position in a power outage
- Weight capacity 450 LBS
- Semi-Electric beds: adjustable head and foot sections, overall height it manually operated with a crank
- Full Electric Beds: all 3 functions (head, foot, overall height) are operated with a push-button pendant

**HOSPITAL BED RAILS**



**HALF RAILS**

- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame
- Attractive brown-vein finish is easy-to-maintain and less institutional looking
- Constructed of 1" steel



**FULL RAILS**

- Attractive brown-vein finish
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove

**ACCESSORY ITEMS**

**HOSPITAL BED FITTED SHEET 36X80'**  
 50-50 COTTON POLYESTER \$23.95 -EA

**THERMAL BLANKET 100% COTTON**  
 66x90" \$25.50 -EA

**HOSPITAL BED FITTED SHEET SET**  
 1-FITTED,1-FLAT,1-PILLOWCASE  
 50-50 COTTON POLYESTER \$49.95 -EA



**BLANKET SUPPORT**  
 FITS ALL BEDS, ADJUSTS TO 23" HIGH  
 ALUM/BRONZE FINISH  
 KEEPS BLANKET OFF FEET \$29.95 -EA

**BED RAIL BUMPER PADS**  
 72Lx14Hx1.5" \$79.95 -PR






**THERAPEUTIC 5-ZONE SUPPORT MATTRESS**  
 Deluxe horizontal, cross-cut foam mattress provides comfort, support, and pressure redistribution over 5 therapeutic pressure zones

- Better compliance and durability than innerspring mattresses allow for improved patient satisfaction and fewer service calls
  - Nonskid bottom keeps mattress safely in place
  - Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633  
 15019, \$329.95

**GRAVITY-7 MULTI-LAYERED MULTI-ZONED MATTRESS**

High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction



- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
  - Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
  - Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
  - High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
  - Visco foam heel section offloads pressure from vulnerable heels
  - Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633  
 15770, \$369.95




**ALTERNATING PRESSURE OVERLAY SYSTEM**

LVA1004 Complete System, \$135.50  
 LVA2030BGE Pad only, \$34.50

**Comfortable Weight Distribution**

Provides relief through alternating inflation and deflation



**MED-AIRE 8" ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM**

Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control.  
 36Wx80Lx8H" Fluid resistant nylon cover



- Designed to prevent, treat and heal pressure ulcers in the home or long-term care setting
- Control unit offers 10-minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low

- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
  - 8" deep air cells are easily removed and replaced with air leakage plugs
  - Removable Nylon stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
  - CPR valve allows for rapid deflation
  - Tubing quickly and easily disconnects with EZ Lock/Quick Release design
- 14027, \$1299.00

**Premium Guard Gel Overlay**

Best in class gel overlay optimizes pressure redistribution and is now available in up to 500lbs weight capacity

- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration



14893 PremGuardGel/FoamOverlay34x76x3.5" 300lb Capacity, \$269.95  
14901 PremGuardGel/FoamOverlay42x76x3.5" 500lb Capacity, \$495.00

<b>HOSPITAL BED PACKAGE PRICING, DELIVERY AND SET-UP NOT INCLUDED*</b>				
<b>PKG#</b>	<b>Bed Type</b>	<b>Rail Type</b>	<b>Mattress</b>	<b>Price</b>
SF5Z	SEMI-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,399.00
SFG7	SEMI-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,450.00
SH5Z	SEMI-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,399.00
SHG7	SEMI-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,450.00
FF5Z	FULL-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,699.95
FFG7	FULL-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,749.95
FH5Z	FULL-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,699.95
FHG7	FULL-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,749.95

