


| Description | Add'l Info | Item # | Price |
|-------------|------------|--------|-------|
|-------------|------------|--------|-------|

PREMIUM ANEROID SPHYGMOMANOMETER WITH MATCHING CASE



Features:
 Latex free, Nylon cuff, Index and Range Markings, Artery indicator mark, Owner ID label, Gauge holder, Standard inflation bulb and Air release valve, Matching 6"x9" nylon carrying case
 Cuff Range: 10" - 16" arm circumference
 Warranty: Five years with lifetime gauge calibration

- 882-ROY.....ROYAL BLUE
- 882-HUN.....HUNTER GREEN
- 882-NAV.....NAVY
- 882-PUR.....PURPLE
- 882-BLK.....BLACK
- 882-BUR.....BURGUNDY
- 882-AQS.....AQUA SEA
- 882-CBL.....CIEL

| | | | |
|-----------------------------------|------------------|---------------|-------|
| ANEROID SPHYG BLACK NYLON CUFF LF | LG ADULT 34-50CM | ADC 720-12XBK | 85.95 |
| ANEROID SPHYG BLACK NYLON CUFF LF | ADULT 23-40CM | ADC 760-11ABK | 48.00 |
| ANEROID SPHYG BLACK NYLON CUFF LF | CHILD 13-19.5CM | ADC 760-9CBK | 48.00 |

A&D AUTOMATIC BLOOD PRESSURE MONITORS



Product Features

- Clinically Validated for Accuracy
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- One Button Operation
- Store 60 Readings each for up to 4 users
- Average reading feature provides trending
- Time & date stamp for every reading
- AccuFit Plus Cuff fits most arms 8.6-16.5"

PREMIUM WRIST BLOOD PRESSURE MONITOR

Advanced, multi-user blood pressure monitor that provides trusted accuracy and reliable measurements every time **UB543, \$79.95**

UA767FAC (WITH AC ADAPTER).....\$95.50

UA767F (WITHOUT AC ADAPTER).....\$85.95

| | | | |
|--|----------------------------|---------------|---------------|
| DIGITAL BP MONITOR AUTO XL CUFF 16.5-23.6" | WITH MEMORY,AC ADAPTER | UA789AC | 189.95 |
| DIGITAL BP MONITOR AUTO SMALL ADULT CUFF | WITH MEMORY,AC ADAPTER | UA767PSAC | 95.50 |
| OTOSCOPE,POCKET XENON 3X MAG | INCLUDES:3-SPECS,2-AAA BAT | PRESTIGE 3470 | 69.95 |
| PENLIGHT WITH PUPIL GAUGE DISP | EACH,6/PK | ADC 351 | 2.95 |
| PENLIGHT REUSE METALITE WHT PUPIL | 2-AAA BAT INCLUDED | ADC 352WP | 6.95 |

3470 XENON POCKET OTOSCOPE





PULSE OXIMETERS

Fingertip measurement of SpO2% (Oxygen Saturation) and pulse rate. Measures quickly and accurately. LED display screen can display information in two directions. Low power consumption with auto shut-off feature. Up to 30 hours of continuous use. Lightweight and compact design.

- 459-BLK BLACK**
 - 459-HKW HOT PINK**
 - 459-PUR PURPLE**
 - 459-REW RED/WHITE**
 - 459-ROW BLUE/WHITE**
- \$36.50**



MODEL-S PULSE OXIMETER
\$29.95

| Description | Add'l Info | Item # | Price | | |
|---|---|-------------------------------------|-------------------------|---------------|--------------|
| STETHOSCOPES | | | | | |
|  | LITTMANN CARDIOLOGY - Diagnose and make decisions with confidence <ul style="list-style-type: none"> • Hear subtle changes in patient status • Useful in critical care and challenging environments such as the ED, ICU, Cardiac ICU, step-down unit and other dynamic locations • Ergonomic and high profile construction; 40% larger chestpiece and 60% deeper bell than the 3M™ Littmann® Classic III™ Stethoscope • Designed for use with adult and pediatric patients • Tunable, dual-sided stainless steel chestpiece with open or closed bell • Non-chill bell sleeve for greater patient comfort • Soft-sealing eartips provide optimal comfort and excellent sound occlusion • Made in the USA without natural rubber latex | | | | |
| | STETHOSCOPE,LITTMANN CARDIOLOGY 4 | NAVY TUBE 27" | 3M 6154 | 279.95 | |
| | STETHOSCOPE,LITTMANN MASTER CARDIOLOGY | BLACK TUBE 27" | 3M 2160 | 339.25 | |
|  | Littmann's latest addition to their line. The Littmann Classic III features patented tunable technology that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece. High acoustic sensitivity for performing general physical assessments. The bell stays clear of dirt and debris by keeping it covered with the small diaphragm. | | | | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | BLACK TUBE 27" | 3M 5620 | 149.95 | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | NAVY TUBE 27" | 3M 5622 | | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | CARIBBEAN BLUE TUBE 27" | 3M 5623 | | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | BURGUNDY TUBE 27" | 3M 5627 | | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | CIEL BLUE TUBE 27" | 3M 5630 | | |
| | STETHOSCOPE LITTMANN CLASSIC 3 | LAVENDER TUBE 27" | 3M 5632 | | |
| |  | STETHOSCOPE LITTMANN LITWEIGHT BLK | BLACK TUBE 28" | 3M 2450 | 72.50 |
| | | STETHOSCOPE LITTMANN LITWEIGHT BRG | BURGUNDY TUBE 28" | 3M 2451 | |
| | | STETHOSCOPE LITTMANN LITWEIGHT CEIL | CEIL BLUE TUBE 28" | 3M 2454 | |
| | | STETHOSCOPE LITTMANN LITWEIGHT CRB | CARRIBEAN BLUE TUBE 28" | 3M 2452 | |
| | | STETHOSCOPE LITTMANN LITWEIGHT PINK | PINK TUBE 28" | 3M 2456 | |
|  | STETHOSCOPE CLINICAL PLUS LTW BLK | 31",4.9oz SS STEM | S123-BLK | 39.95 | |
| | STETHOSCOPE CLINICAL PLUS LTW PINK | 31",4.9oz SS STEM | S123-HPK | | |
| | STETHOSCOPE CLINICAL PLUS LTW NAVY | 31",4.9oz SS STEM | S123-NAV | | |
| | STETHOSCOPE CLINICAL PLUS LTW PURPLE | 31",4.9oz SS STEM | S123-PUR | | |

The Adtemp™ 418N features:

- 8-second measurement using proprietary predictive technology
- Range 93.2°F to 111.0°F ±.2°F or 34.0°C to 43.9°C ±.1°C depending upon scale selection
- Large LCD display with color-coded backlighting (green for readings in the normal range and red for readings above 99.5°F/37.5°C)
- Replaceable 3V battery (CR2032) provides up to 2,000 measurements
- Measures 5.31"L x 1.14"W x 0.59"H, weighs 0.77 oz (22 g) with battery
- Integral scabbard-type carry case, Includes five probe sheaths
- Meets the accuracy requirements for intermittent thermometers under ASTM E1112

ADC 418N, \$16.95



NON-CONTACT THERMOMETERS



- Non-Contact, Fast , Accurate Readings
 - Measures Body, Surface and Ambient Temperatures
 - Backlit Digital LCD Display Screen
 - High Temperature Alarm
 - Measures in Both Fahrenheit & Celsius
 - Features 32 data memory
- Includes two AAA batteries, 1 year warranty

DT-29, \$49.95