

 <p>VERSATILT OVERBED TABLE 15x30" WALNUT MELANINE TOP 26.5-39" INFINITE HT ADJUST 950H, \$249.95</p>	 <p>OVERBED TABLE WITH LOCKING WHEEL STOPS NON-TILT TOP:15.5Wx30L LVA1022, \$149.95</p>	 <p>Omni Tray with Bamboo 20x15" Swivel Top STAN 2600, \$229.95</p>	 <p>SILVER VEIN FINISH 300LB CAP FOOT STOOL, \$44.50 FOOT STOOL+HANDLE, \$65.50</p>
---	---	---	---

PATIENT LIFTS



- Silver Vein steel construction
- 6-point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base
- 450 LB weight capacity
- Lifetime warranty (frame), 1 year limited (pump)
13023SV, \$895.00



The PL400HE is Bestcare’s answer to home care needs. Featuring a 400 LB maximum weight capacity and Bestcare battery powered control system for power lifting, the PL400HE is easily adaptable for a wide range of challenges while still maintaining a high standard of comfort for patients. The use of Bestcare’s battery powered electronic system provides a low-cost solution for power lifting while providing convenience for caregivers. Designed for home care use and effective at providing quality care.

- 400 LB carrying capacity
- Six-point spreader bar
- Electronic control system
- Optional digital weigh scale
- Optional two-point spreader bar
PL400HE, \$1,975.00

SLINGS



MESH SLING FULL BODY

WITH COMMODE OPEN

SL-R110	MED	WEIGHT RANGE 100-195 LBS 450LB CAP
SL-R111	LG	WEIGHT RANGE 150-300 LBS 450LB CAP
WITH COMMODE OPENING:		
SL-R114	MED	WEIGHT RANGE 100-195 LBS 450LB CAP
SL-R115	LG	WEIGHT RANGE 150-300 LBS 450LB CAP
SL-R116	XLG	WEIGHT RANGE 210-400 LBS 450LB CAP

SLING DIVIDED LEG MESH

HOYER TYPE

SL-UM882	MED	FITS WEIGHT RANGE 140-200 LBS
SL-UM883	LG	FITS WEIGHT RANGE 200-450 LBS
HOYER TYPE:		
SL-HC70002	MED	U-SLING 99-210 LBS 500LB CAP
SL-HC70001	LG	U-SLING 150-300 LBS 500LB CAP
SL-HC70000	XLG	U-SLING 270-600 LBS 600LB CAP
SL-HC70012	MED	W/HEAD SUP'RT 99-210 500LB CAP



FREEDOM

TRANSFER PATIENT LIFT



ERGONOMIC HANDLES

Comfortable and secure hand grips in the front and rear for easy transport.

01

SEAT & FRAME LOCK

Dual locking to ensure the safety of the patient and the caregiver.

02

DUAL CONTROLLED

Features built-in control and a detachable remote for lifting and lowering the patient.

03

ADJUSTABLE

The seat height may be adjusted from 15.5" to 23.5" to accommodate height differences when transporting.

04

FOOTREST

Keep the patients feet comfortable and secure when transporting patients.

05



THE
BED

THE
TOILET

THE DINING
ROOM

THE LIVING
ROOM

YOUR
FAMILY

\$2,799.00

Overall Height: 34.5" – 42.5"
 Seat Width: 19.5"
 Seat Height: 15.5" – 23.5"
 Commode Pail: 6 Qt
 Required Clearance: 4.5"
 Rear Wheel Size: 3"
 Wheel Base Length: 25"
 Wireless Lithium Battery
 Footrest Size: 13.75" x 6"
 Wheel Base Width: 22.3"



TRAPEZE BAR

Frame is constructed of heavy-gauge, chrome plated, steel tubing for maximum strength. Wall bumpers and padded mounting brackets protect bed ends and wall surfaces. Assembles without tools; mounts to bed. Optional base (shown) allows trapeze to be free standing. 250 lb (113 kg) assisted weight capacity.

MDS80615T (TRAPEZE ONLY), \$195.95

MDS80616T, (FLOOR BASE), \$169.95

ELECTRIC HOSPITAL BEDS



- Constructed from durable steel with a spring deck and split-frame design
- Spring deck is adjustable from 15.7"-22.8"
- Motor is compatible with a 9-volt battery to adjust bed position in a power outage
- Weight capacity 450 LBS
- Semi-Electric beds: adjustable head and foot sections, overall height it manually operated with a crank
- Full Electric Beds: all 3 functions (head, foot, overall height) are operated with a push-button pendant

HOSPITAL BED RAILS



HALF RAILS

- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame
- Attractive brown-vein finish is easy-to-maintain and less institutional looking
- Constructed of 1" steel



FULL RAILS

- Attractive brown-vein finish
- Constructed of 1" steel
- Telescopes with movement of bed
- Easily adjusts up or down with spring-loaded release
- Spring-loaded crossbars are easy-to-install and remove

ACCESSORY ITEMS

HOSPITAL BED FITTED SHEET 36X80'
50-50 COTTON POLYESTER \$23.95 -EA

THERMAL BLANKET 100% COTTON
66x90" \$25.50 -EA

HOSPITAL BED FITTED SHEET SET
1-FITTED,1-FLAT,1-PILLOWCASE
50-50 COTTON POLYESTER \$49.95 -EA



BLANKET SUPPORT
FITS ALL BEDS, ADJUSTS TO 23" HIGH
ALUM/BRONZE FINISH
KEEPS BLANKET OFF FEET

\$29.95 -EA

BED RAIL BUMPER PADS
36Lx14Hx1.5", \$69.95 -PR
48Lx14Hx1.5", \$95.50 -PR
72Lx14Hx1.5", \$95.50 -PR






THERAPEUTIC 5-ZONE SUPPORT MATTRESS
 Deluxe horizontal, cross-cut foam mattress provides comfort, support, and pressure redistribution over 5 therapeutic pressure zones

- Better compliance and durability than innerspring mattresses allow for improved patient satisfaction and fewer service calls
- Nonskid bottom keeps mattress safely in place
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
 Meets flammability standard CFR 1632 and CFR 1633 15019, \$329.95

GRAVITY-7 MULTI-LAYERED MULTI-ZONED MATTRESS

High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction



- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
 - Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
 - Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
 - High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
 - Visco foam heel section offloads pressure from vulnerable heels
 - Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture. Meets flammability standard CFR 1632 and CFR 1633 15770, \$399.95




ALTERNATING PRESSURE OVERLAY SYSTEM

LVA1004 Complete System, \$135.50
 LVA2030BGE Pad only, \$34.50

Comfortable Weight Distribution

Provides relief through alternating inflation and deflation



MED-AIRE 8" ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM

Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control.
 36Wx80Lx8H" Fluid resistant nylon cover



- Designed to prevent, treat and heal pressure ulcers in the home or long-term care setting
- Control unit offers 10-minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low

- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
 - 8" deep air cells are easily removed and replaced with air leakage plugs
 - Removable Nylon stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
 - CPR valve allows for rapid deflation
 - Tubing quickly and easily disconnects with EZ Lock/Quick Release design
- 14027, \$1299.00

Premium Guard Gel Overlay

Best in class gel overlay optimizes pressure redistribution and is now available in up to 500lbs weight capacity

- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration



14893 PremGuardGel/FoamOverlay34x76x3.5" 300lb Capacity, \$269.95
14901 PremGuardGel/FoamOverlay42x76x3.5" 500lb Capacity, \$495.00

HOSPITAL BED PACKAGE PRICING, DELIVERY AND SET-UP NOT INCLUDED*				
PKG#	Bed Type	Rail Type	Mattress	Price
SF5Z	SEMI-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,399.00
SFG7	SEMI-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,450.00
SH5Z	SEMI-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,399.00
SHG7	SEMI-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,450.00
FF5Z	FULL-ELECTRIC	FULL RAIL	15019 MATTRESS	\$1,699.95
FFG7	FULL-ELECTRIC	FULL RAIL	15770 MATTRESS	\$1,749.95
FH5Z	FULL-ELECTRIC	HALF RAIL	15019 MATTRESS	\$1,699.95
FHG7	FULL-ELECTRIC	HALF RAIL	15770 MATTRESS	\$1,749.95

